



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 8995



FreeZone Stopping Tray Dryer 115V

DESCRIPTION: This stopping tray mounts on top of a 6-, 12- or 18-Liter Console Freeze Dryer (not included). It is used in a laboratory setting to facilitate the lyophilization process. It helps in freeze drying moderate to large loads of vials or serum bottles and stoppering them under vacuum.

FEATURES:

- Three 16.5" × 12.0" adjustable processing shelves. The spacing between the shelves accommodates containers with stoppers up to a maximum height of 75 mm and minimum height of 38 mm when all three shelves are used and increases accordingly when two shelves or one shelf is used
- Three probes for monitoring samples or shelf temperature
- Integral support stand
- Data logging capabilities that can be exported via USB or Ethernet from the FreeZone Console Freeze Dryer
- Stopping mechanism that inflates a diaphragm via atmospheric pressure. No compressed gas required

TECHNICAL SPECIFICATIONS:

Weight (lbs/kg)	423/192
Dimensions, W x D x H (in/cm)	32.4 x 24.6 x 27.1/82.3 x 62.5 x 68.8
Electrical	115V, 60Hz, 16 amps
Conformance	CAN/CSA C22.2 UL
Vacuum Pump Displacement Required	163 L/min or larger displacement

APPLICABLE PRODUCTS:

- 8994: Vacuum Pump, 195L/Min 115V (Freezone Sys) **(REQUIRED)**
- 8996: FreeZone Plus Freeze Dry Sys. 6L Cascade **(REQUIRED)**

WARRANTY: 1 year

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM