



AirClean® PowderSafe™ Bulk Powder Handling Enclosure, 5'

With Base Stand

Description: When effective powder containment is needed during bulk manipulation, the PowderSafe™ Bulk Powder Handling Enclosure is an economical and space-saving alternative to larger walk-in enclosures. With a unique drum access port incorporated into the base of the enclosure, the drum can be raised and sealed to the access port, effectively extending the enclosure's containment area into the drum. Designed to provide a turbulence-free airflow environment, PowderSafe Bulk Powder Handling enclosure moves air in a horizontal pattern to maximize containment while minimizing sample loss and balance instability.

Features:

- Dark blue base for ease of cleaning
- Disposal port for weighing vessels
- HEPASafe™ filter change-out
- AirSafe™ automatic safety controller
- Bulk container access port cover
- Vapor-proof fluorescent light
- No installation required – No ductwork required
- White epoxy coated stand with leveling feet
- Electrical and utility cord access ports
- Optional hydraulic lift for bulk powder drum

Technical Specifications:

Weight (kg/lb)	226.8 / 500.0
-----------------------	---------------

USA | Canada | Australia | International
medisca.com

External Dimensions : W x D x H (cm/in) *height includes base stand	152.4 x 91.4 x 236.2 / 60.0 x 36.0 x 93.0
Electrical	110 V, 4.3 A
Face Velocity (lfm)	40 - 120
Airflow Monitor	AirSafe™ Automatic Safety Controller
Blower Requirements	Brushless-Sparkless
Lighting	Vapor-Proof Fluorescent
Sash Height (cm/in)	20.3 / 8.0
Internal Height (cm/in)	88.9 / 35.0
Conformance	UL, CSA, ASHRAE 110 Tested, ANSI Z9.5

Associated Products:

- Pre-filter: 7278 – AirClean® PowderSafe™ Type A and Bulk Pre-Filter, 3’/5’ (4 x Required)
- Main Filter(s): 7277 – AirClean® PowderSafe™ Type A, B, and Bulk HEPA Filter, 3’/5’ (4 x Required)
- Drum Lift: 6946 – AirClean® PowderSafe™ Bulk Powder Enclosure Manual Drum Lift

Warranty: 1 year

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM

Please note that the information contained herein is for informational purposes only. The information is based on recognized compendia, manufacturers' instructions and/or other authoritative sources and texts. MEDISCA takes no responsibility for the validity or accuracy of this information. Customer MUST refer to USP, FDA, Health Canada, TGA and/or their State Board to ensure they comply with the specific regulations and ensure that appropriate procedures are in place.