

SINCE
2015

MD[®] DISPENSERS

A Brand Built on Elegance, Accuracy and Versatility



 MEDISCA[®]

The Next Generation of MD[®] Dispensers

MEDISCA has expanded and diversified the MD Dispensers portfolio to cater to the needs of different niche markets. This line provides improved efficiencies, design differentiation and functional benefits that ultimately enhance your patients' overall product experience.

Get The MD Facts



First of its kind, diverse airless dispensing portfolio with the convenient UV-resistant view window



Portfolio has expanded from 30 to 115 dispenser combinations, based on customer feedback



Used by Over 50,000 patients worldwide



Withstand Air Pressure, which makes it ideal for travelling



UV-resistant view window proven to block 99% of UVA and UVB rays

VERSATILITY REDEFINED

MD® JARS

When Elegance Meets Accuracy

MD Jars feature a refined design and cosmetic appeal. With a convenient “press to use” mechanism, MD Jars allow for hygienic use and protection from oxidization without sacrificing accurate measuring of the dose. Available in 50 mL size with 0.5 mL actuation.

Suggested Applications:



Cosmetic Applications



Skin Conditions



Get Your Sample Formula

F 004 625v2: Hydroquinone 2%,
Tretinoin 0.015% Topical Cream



MEDISCA
NETWORK

www.medisca.net

For further product details and sizing please see the reference chart on page 11.

MD® PUMPS

High Performance in a Pump

MD Pumps are engineered keeping accuracy and precision at the forefront without compromising pharmaceutical elegance. Featuring a UV-resistant view window, this dispenser ensures that patients can easily and accurately dose medications. Available in a range of sizes: 15 mL to 240 mL.

Suggested Applications:



HRT
Applications



Pain Management
Applications



Skin
Conditions



Sunscreen/Cosmetic
Applications

For further product details and sizing please see the reference chart on page 10.



Voice of Customer

"MD Dispensers provide an ease-of-use solution for applications. A see-through indicator helps patients know when their medication is low so it can be ordered without a break in treatment. MD Dispensers are the only device that topical preparations should be dispensed in."

- Ron Yochim, Compounding Pharmacy Manager
Cook's pharmacy, Canada



NOTE: MD Actuators are sold separately from MD Pumps.

MD[®] ACTUATORS

Precision Right at Your Fingertips

MD Actuators' exterior design features various colors to facilitate medication identification. The Actuators fit all MD Pump bottles for added flexibility of metered dosing options. Available in green, blue, pink and white and 4 dosage options: 0.25mL, 0.5 mL, 1 mL and 1.5 mL.

Suggested Applications:



HRT
Applications



Pain Management
Applications



Skin
Conditions



Sunscreen/Cosmetic
Applications




**On-Demand
Webinar**

Compounding Opportunities in
Transdermal Delivery



www.lp3network.com



For further product details and sizing please see the reference chart on page 10.

MD[®] SYRINGES

Not Your Ordinary Syringe, Even for Little Troubles

With a sleek medical look, a silicone tip for smooth application and excellent evacuation, the MD Syringes are elegantly designed and easy to operate. They provide a repeatable and accurate dose each time for optimal efficiency and targeted delivery. Available in 5 mL and 10 mL size with 0.15 mL actuation.

Suggested Applications:



Localized
Scarring



Skin
Conditions



Veterinary
Applications



Cosmetic
Applications



On-Demand Webinar

Optimizing the Success of Wound
Healing and Scar Reduction



www.lp3network.com



Silicone tip for
smooth application

Medical look

See-through
UV-resistant
view window

Close/lock
feature allows
for easy travel

For further product details and sizing please see
the reference chart on page 11.

MD® PENS

A Vet's Best Friend

The MD Pen features a slim pen-type look that facilitates the application of the medication. Unlike a felt tip, a silicone tip allows for smooth application eliminating the risk of bacterial build-up or contamination. Available in 6.5 mL size with 0.10 mL or 0.15 mL actuation.

Suggested Applications:



Localized Scarring



Skin Conditions




Under Eye Cream/Acne



Insect Bites



Veterinary Applications



Get Your Sample Formula

F 006 918: Methimazole 5mg/0.1 mL Transdermal PLO Gel



www.medisca.net

For further product details and sizing please see the reference chart on page 11.





Voice of Customer

"We recently launched a compounding business at our pharmacy and MEDISCA has been instrumental in getting us to where we are today! We especially love MEDISCA's Foamil™ Base. It is easy to work with and the MD Foamers dispense a silky foam that our customers love."

- Sara Kemp, BSc Pharm, RPh
HealthConnect Pharmacy, Canada

MD® FOAMERS

Get a Silky Foam That Customers Will Love

Did you know that the properties of the foam produced are dependent on both the quality of the foam base and the dispenser used? The MD Foamer is essential to creating a dense and creamy foam. Available in 50 mL size with 0.5 mL actuation or 110 mL with 0.7 mL actuation.

Suggested Applications:



Face Wash
or Cleanser



Hair Care
Treatments



MD[®] MINIS

Pump and Go!

MD Minis are compact, practical and are specifically developed for patients with on-the-go busy lifestyles, allowing them to discreetly apply their medications. Available in 5 mL and 10 mL sizes with 0.15 mL actuation.

Suggested Applications:



HRT
Applications



Pain Management
Applications



Sampling to
Physicians




Sunscreen/Cosmetic
Applications

Colored bands allow for easy
medication identification

Airless dispensing

Metered dosing

See-through
UV-resistant
view window



Get Your Sample Formula

F 007 841: Progesterone 25 mg/0.15 mL
Topical Cream

F 007 842: Baclofen 2%, Diclofenac
Sodium 3% Topical Cream



MEDISCA
NETWORK

www.medisca.net

For further product details and sizing please
see the reference chart on page 11.



MD® DISPENSERS REFERENCE CHART

MD Pumps*

PRODUCT	PROD. NO.	SIZE
	6337	15 mL
	7371	30 mL
	7372	50 mL
	7373	80 mL
	7374	100 mL
	7375	150 mL
	7376	200 mL
	7216	240 mL



MD Primer

PRODUCT NO.: 6017

Manual priming device that will be used for MD pumps. The MD Primer will push air in from the bottom of the pump in order to raise the pumps internal piston.

MD CapSnapper

PRODUCT NO.: 7085

Once MD Dispensers are filled up to the designated volume, snapping is quick and effortless with the use of the CapSnapper.




MD Actuators*

PRODUCT	PROD. NO.	METERED DOSE
	A 6344	0.5 mL
	B 6341	1.0 mL
	C 6338	1.5 mL
	D 5947	0.25 mL
	A 6346	0.5 mL
	B 6343	1.0 mL
	C 6340	1.5 mL
	D 5945	0.25 mL
	A 6345	0.5 mL
	B 6342	1.0 mL
	C 6339	1.5 mL
	D 5946	0.25 mL
	A 7370	0.5 mL
	B 7369	1.0 mL
	C 7217	1.5 mL
	D 5944	0.25 mL
	S 6688	1.0 mL

*MD Actuators and Pumps are sold separately.

MD® DISPENSERS REFERENCE CHART

MD Minis

PRODUCT	PROD. NO.	SIZE	METERED DOSE
	6252	5 mL	0.15 mL
	6248	10 mL	0.15 mL
	6254	5 mL	0.15 mL
	6250	10 mL	0.15 mL
	6253	5 mL	0.15 mL
	6249	10 mL	0.15 mL
	6255	5 mL	0.15 mL
	6251	10 mL	0.15 mL

MD Dispensers

PRODUCT	PROD. NO.	SIZE	METERED DOSE
	7080	6.5 mL	0.10 mL
	7377	6.5 mL	0.15 mL
	7379	5 mL	0.15 mL
	7378	10 mL	0.15 mL
	6701	50 mL	0.5 mL
	6700	110 mL	0.7 mL
	7380	50 mL	0.5 mL

Discover The Right Dispensing Solution



PERSONALIZED MEDICATION
Starts With the Perfect Partnership.



Products
www.medisca.com



Services
www.medisca.net



Education
www.lp3network.com