



Medisca Transfer Pipette

Description: The Medisca Transfer Pipette is ideal for preparation of NovaFilm™ Oral Transmucosal Films as well as other applications that require repetitive and precise dispensing. Comes equipped with two 50 mL syringe tips and one adapter to connect the tip to the Medisca Transfer Pipette.

Features:

- Adjusts to multiple dispensing volumes – from 1 mL to 5mL using the 50 mL tip
- Facilitates the transfer of liquids with ranging viscosity and density
- Can deliver up to 48 dispensing steps in succession
- Tips are made from high quality materials without soft lubricants or other additives
- Continuously adjustable
- Slim, light and robust design
- Easy to use

*Note: the first and last dispensing liquid should be discarded for precision.

Technical Specifications:

Material	
Piston	HDPE
Cylinder	PP

USA | Canada | Australia | International
medisca.com

Dispensing Volume:

Dial Setting		1	2	3	4	5		
Number of Steps		48	23	15	11	8		
Product Code	Tip Volume	Sample Volume					Accuracy error	Precision error
5497	50.mL	1000 µL (1 mL)	2000 µL (2 mL)	3000 µL (3 mL)	4000 µL (4 mL)	5000 µL (5 mL)	± 0.2%	≤ 0.2%

Applicable Products:

- NovaFilm™ Gel Base ([3187](#))
- Blister Mold, NovaFilm™ OTF ([5524](#))
- Rack for NovaFilm™ OTF Blister Mold ([5523](#))
- Syringe Tips, 50 mL, w/1pc Adaptor, For Medisca Transfer Pipette ([5497](#))