



**Mixer, Mazerustar KK-400W, Counter Balance Weight, 25 g**

**Description:** This counterbalance weight can be placed directly in the MAZ<sup>®</sup> KK-400W holder facilitating the counterbalance step required before initiating any cycle on the MAZ<sup>®</sup>. The corresponding mass engraved on each weight facilitates identification.

**Technical Specifications:**

Material	Acetal (Sustarin C)
Color	Natural (White)
Outer Diameter (mm)	70
Inner Diameter (mm)	18
Height (mm)	9
Mass (g)	24 - 26

USA | Canada | Australia | International  
medisca.com

Please note that the information contained herein is for informational purposes only. The information is based on recognized compendia, manufacturers' instructions and/or other authoritative sources and texts. MEDISCA takes no responsibility for the validity or accuracy of this information. Customer MUST refer to USP, FDA, Health Canada, TGA and/or their State Board to ensure they comply with the specific regulations and ensure that appropriate procedures are in place.