



Ohaus Dot Matrix Printer SF40A, with Cable

Description: Introducing the SF40A, a portable printer that combines high performance and versatility with Ohaus's trusted quality and durability. This user-friendly printer offers advanced functionality at an exceptional value. Ideal for use alongside Ohaus balances and scales, the SF40A is perfect for laboratory, industrial, and educational applications, ensuring accurate documentation of weighing results. Enhance your workflow with this reliable tool designed for precision and ease of use.

Features:

- Paper roll, ink ribbon, cable, cable adapter, and power cable included.
- The statistics and totalization modes and real time clock complement Ohaus weighing instruments' application modes and offer functionality that assists the reporting processes.
- SF40A helps support the efficiency of the weighing and recording processes. RS232 interface and PC configuration software make setting up and using the printer fast and simple.
- With the ability to print in 13 different languages, the SF40A can print in the language of the connected device to support universal operation: English, Spanish, German, Italian, French, Chinese, Korean, Japanese, Russian, Polish, Portuguese, Hungarian, and Czech.

Technical Specifications:

| | |
|--|-------------------------------------|
| Printer dimensions: w x d x h (cm/in) | 14.0 x 20.5 x 7.5 / 5.5 x 8.1 x 2.9 |
| Net weight, paper roll included (kg/lb) | 0.8 / 1.8 |

USA | Canada | Australia | International
medisca.com

| | |
|--------------------------------------|--------------------------|
| Interface | RS232 |
| Power requirements | AC adapter (included) |
| Paper roll: w x Ø (cm/in) max | 5.75 x 6.0 / 2.26 x 2.36 |
| Print speed | 1.2 lines per second |
| Number of columns | 40 |

Associated Products:

- 6238 – Ohaus Printer Paper Roll, SF40A
- 7186 – Ohaus Ink Ribbon Cartridge, SF40A

Warranty: 1 year

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM

Please note that the information contained herein is for informational purposes only. The information is based on recognized compendia, manufacturers' instructions and/or other authoritative sources and texts. MEDISCA takes no responsibility for the validity or accuracy of this information. Customer MUST refer to USP, FDA, Health Canada, TGA and/or their State Board to ensure they comply with the specific regulations and ensure that appropriate procedures are in place.