



## **Capsule Machine, ProFiller 3600, Size 000**

**Description:** The ProFiller 3600 is the only Inox (316 stainless steel) hand-held 300 hole filler that is dishwasher safe & can be used under a powder hood or in a glove box. The Change Parts and Orienter can also be used with the ProFiller 3700 bench-top filler for maximum flexibility. Included Parts: Hand-Held Filler Base, Orienter Base, Locking Plate, 60-Pin Tamper, Powder Spreader, Powder Tray, Powder Tray Clamps, Plastic Storage Box, Capsule Bins, Measuring Bowls and Change Parts (Caps Tray, Body Sheet Set (A, B, C), Spacer Plate, Capsule Locker and Orienter Tray).

**\*\*Please note that 000 change parts are 180 holes only. Due to the size of these capsules, the machines are not able to accommodate 300 capsules.**

### **Features:**

- Shipped fully assembled. Ready for use in less than 1 hour
- Easier to use than other 200 or 300-hole systems
- All parts are dishwasher safe. Can be dried in a hot air oven or with a portable hair dryer
- Change capsule sizes in only 5 minutes
- Orienter Base accommodates capsule sizes 00 – 4
- All metal parts Inox (SS316) for reduced static
- Unique features include lightweight 2-piece Orienter, Orienter Locator and Oversize Powder Tray & Push Plate

USA | Canada | Australia | International

[medisca.com](http://medisca.com)

- Designed for fast disassembly & visual verification of cleaning
- Fill weights, USP/EP method within  $\pm$  1-3 % of average
- Works with all capsules including gelatin, vegetarian/HPMC and starch capsules
- Tamper size range optimized for higher fill weights and faster filling
- Indicator gauges to assist in capsule orientation, separation and locking
- Fits under most powder hoods and in glove boxes
- 3600 fully locks capsules after filling, unlike most hand-held fillers
- Compliance: CE, cGMP

**Warranty: 1 year**

USA | CANADA | AUSTRALIA PACIFIC | EUROPE  
WWW.MEDISCA.COM

Please note that the information contained herein is for informational purposes only. The information is based on recognized compendia, manufacturers' instructions and/or other authoritative sources and texts. MEDISCA takes no responsibility for the validity or accuracy of this information. Customer MUST refer to USP, FDA, Health Canada, TGA and/or their State Board to ensure they comply with the specific regulations and ensure that appropriate procedures are in place.