



## BD PhaSeal™ Connector Luer Lock C45 for N35 Connection (Sterile)

**Description:** PhaSeal™ is a system for the safe handling of hazardous drugs. As a clinically proven closed-system drug transfer device (CSTD), it has been proven to prevent hazardous drug exposure from preparation and administration to waste disposal.

The Luer Lock Connector C45 attaches to the patient's IV line, providing a closed connection of hazardous drugs to the Injector. This connector is designed for use with needle-free ports with internal pistons. The C45 mates with the BD PhaSeal Injector Luer Lock (N35) and provides a route for safe drug transfer via injection or infusion.

### Features:

- The only CSTD with extensive third party evidence proving:
  - Prevention of drug exposure
  - Reduction of surface contamination
  - Elimination of human uptake
  - Cutting drug waste to enable savings
- Protects the microbial state of vial contents\*:
  - May reduce impact of drug shortages
  - Helps fund investments in healthcare worker safety
  - May enable better healthcare outcomes

*\*Data on file at BD*

USA | Canada | Australia | International

medisca.com

- Setting the standard in hazardous drug safe handling:  
Validated by 3rd party studies  
Utilized by top NCCN and NCI facilities
- Leading the way:  
First to obtain ONB clearance  
First with clearance to prevent microbial ingress  
First with indications aligned to the NIOSH definition of a CSTD
- Disposable
- BPA-Free
- Latex-Free
- Sterilized using Ethylene Oxide
- Priming Volume (mL): 0.1

**Material Specifications:**

<b>Housing</b>	<b>PP</b>
<b>Membrane</b>	<b>TPE</b>
<b>Protective Cap</b>	<b>PE</b>

**Applicable Products:**

- **7105: BD PhaSeal™ Injector, N35**

USA | CANADA | AUSTRALIA PACIFIC | EUROPE  
WWW.MEDISCA.COM

Please note that the information contained herein is for informational purposes only. The information is based on recognized compendia, manufacturers' instructions and/or other authoritative sources and texts. MEDISCA takes no responsibility for the validity or accuracy of this information. Customer MUST refer to USP, FDA, Health Canada, TGA and/or their State Board to ensure they comply with the specific regulations and ensure that appropriate procedures are in place.