



Syringe Tips, 50mL, w/1pc Adaptor, For Medisca Transfer Pipette

Description: 50 mL syringe tips that go with the Medisca Transfer Pipette. Each syringe tip comes equipped with 1 pc Adaptor in order to connect the tip to the Medisca Transfer Pipette.

Features:

- Non-sterile
- Tips are made from high quality materials without soft lubricants or other additives

Dispensing Volume:

Dial Setting		1	2	3	4	5		
Number of Steps		48	23	15	11	8		
Product Code	Tip Volume	Sample Volume					Accuracy error	Precision error
5497	50.mL	1000 µL (1 mL)	2000 µL (2 mL)	3000 µL (3 mL)	4000 µL (4 mL)	5000 µL (5 mL)	± 0.2%	≤ 0.2%

Applicable products:

- Medisca Transfer Pipette (5531)

USA | Canada | Australia | International
medisca.com

Please note that the information contained herein is for informational purposes only. The information is based on recognized compendia, manufacturers' instructions and/or other authoritative sources and texts. MEDISCA takes no responsibility for the validity or accuracy of this information. Customer MUST refer to USP, FDA, Health Canada, TGA and/or their State Board to ensure they comply with the specific regulations and ensure that appropriate procedures are in place.