



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

**PRODUCT NO.: 6514**






## **CHEMO PREP MAT, COVIDIEN**



16" x 22"

**DESCRIPTION:** This prep mat provides additional safety when compounding or administering chemotherapy drugs. It helps to absorb small spills which may contaminate the work surface and also to prevent drug contact with the patient's skin.

### **FEATURES:**

-  Latex-free: No additional exposure to latex to help reduce possibility of latex allergies
-  Highly absorbent design: Economical method to help control spilled fluids
-  Non-linting: Helps to limit particulates in laminar flow hoods and barrier isolators

### **ADDITIONAL INFORMATION:**

-  It is important to keep the prep mat flat at all times, and away from the front and rear grilles of the laminar flow hood or barrier isolator, to prevent interference with the airflow
-  Because of possible contamination, according to ASHP guidelines, prep mats must be replaced and discarded after preparation of each batch and frequently during extended batch compounding.