

OleaBase™ Plasticized Light Anhydrous and Occlusive

Medisca OleaBase Plasticized Light is a soft, anhydrous gel base of medium viscosity that is suitable for both oral and topical applications and compatible with lipophilic and hydrophilic ingredients. An extension of the original OleaBase Plasticized, OleaBase Plasticized Light is of lower viscosity, is less compatible with salts, and is suitable for oral ingestion.



Backed by Science

- Tested water activity of 0.203 (180 day BUD).
- Suitable for topical and oral use in humans and animals.
- Compatible with lipophilic and hydrophilic ingredients.



Learn More
Product No.: 3286

Features & Benefits

Appearance

Translucent, shiny and elegant oily gel with medium viscosity

Ingredients

- Free of parabens, petrolatum, alcohol, fragrances, and lanolin
 - Contains mineral oil
-

Suitability

- Suitable for use in humans and animals.
 - Suitable for oral and topical applications, including vaginal and rectal mucosal membranes.
-

Compatibility

Compatible with both lipophilic and hydrophilic actives.

Physical Characteristics

Physically stable at 40 °C

Water Activity

0.203

Beyond Use Date

Anhydrous dosage form with a default beyond-use date of 180 days according to USP guidelines.

Formulation Services – Core Support

- 9,000+ Formulas
- 360 Compounding Support Minutes / Year
- 5 Customized Formula Requests / Year

[Access Now](#)

LP3 Network Non-Sterile Training



Gain practical skills in preparing common dosage forms. [Learn More](#)

For More Information

1-800-932-1039 medisca.com

Suggested Products

Medisca MAZ® Mixers [Learn More](#)

Versatile planetary mixer for consistent and efficient mixing, milling, deaeration, and melting needs.



Samix® [Learn More](#)

Versatile and reliable electronic mortar and pestle for your semi-solid preparations.



Medisca MD® Dispensers [Learn More](#)

Built on elegance, accuracy, and versatility.

