



Autoclave Accessories, Sterilization Tape

Description: This autoclave indicator tape has green stripes that turn black during steam sterilization indicating that the pack has been processed. The synthetic, solvent-free adhesive provides secure bonding, yet provides easy removal without traces of residue. High contrast color change ink shows a clear, readily interpreted result, without requiring any guesswork.

Features:

- Black stripes appear after exposure to steam sterilization conditions
- Sticks securely under high-temperature and high-pressure conditions
- Removes easily without leaving sticky residue behind
- Can be used in temperatures ranging from -23°C to +121°C
- Lead-free material
- 54.9m (180ft)
- ISO 11140

USA | Canada | Australia | International

medisca.com