



PRODUCT NO.: 6840



Applicator, Topi-Click, 35ml, 0.5ml per 2 clicks

3 Port, Pink

Description: The Topi-CLICK 3Port uses the same technology as the Topi-CLICK – but with extra ports for faster and easier CLICK-Dosing of thicker medications. The advantages for the patient are added convenience and ease-of-use when dealing with the thicker and refrigerated medications. The NEW Topi-CLICK 3Port is a smart, effective topical applicator that was designed to overcome the various limitations associated with tubes, tubs, syringes and jars.

FEATURES:

- Each click moves the plunger approximately 0.25ml of volume space, thereby dispensing 0.25ml of cream on average
- The most common number of clicks used for one dose is 2 clicks that move the plunger 0.5ml volume space
- Filling with approximately 35ml the Topi-CLICK will provide approximately 64 doses
- Around 0.75ml is used for initial priming
- Can be labeled like a prescription vial – simplified labeling for lower cost – and easier to read
- Domed applicator pad for application without dispensing on to hands first
- Audible, visual and tactile click registration for user recognition
- A “Refill” reminder system helps patients to reorder in a timelier manner
- Designed to resist accidental discharge in purse, pocket or luggage
- BPA and Phthalate free | UV block
- One time use only

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM

TECHNICAL SPECIFICATIONS:

- Dispenser Volume: Up to 37 mL
- Delivery Volume: 5 mL - 30 mL
- Dispense/CLICK: 0.25 mL
- Approx. Residual: 2.5 mL
- Approx. Priming: 0.25 mL or 1 CLICK

USA | CANADA | AUSTRALIA PACIFIC | EUROPE
WWW.MEDISCA.COM