



PRODUCT NO.: 6733



Topi-Click PERL, White, Vaginal Dosing System w/4mL Vaginal Applicator, (35mL, 0.25mL/Click)

Description: The new elegantly designed and comfortable Topi-CLICK Perl is the revolutionary first vaginal dosing system and new vaginal applicator design that is comfortable and accurate, yet easy to use. The Perl was designed with the patients, pharmacist and the physicians in mind.

FEATURES:

- Perl Snap-On Vaginal Applicator is made with clear, flexible soft medical grade polymer
- Automated Topi-CLICK dosing for accuracy and consistency through entire regimen
- Less cost, less clean-up, less waste
- Convenient, fast and easy to use
- Each click moves the plunger approximately 0.25ml of volume space, thereby dispensing 0.25ml of cream on average
- The most common number of clicks used for one dose is 2 clicks that move the plunger 0.5ml volume space
- Filling with approximately 35ml the Topi-CLICK will provide approximately 64 doses
- Around 0.75ml is used for initial priming
- Can be labeled like a prescription vial – simplified labeling for lower cost – and easier to read
- Domed applicator pad for application without dispensing on to hands first
- Audible, visual and tactile click registration for user recognition
- A “Refill” reminder system helps patients to reorder in a timelier manner

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM

- ☑ Designed to resist accidental discharge in purse, pocket or luggage
- ☑ BPA and Phthalate free | UV block
- ☑ One time use only

TECHNICAL SPECIFICATIONS:

- ☑ Dispenser Volume: Up to 37 mL
- ☑ Delivery Volume: 5 mL - 30 mL
- ☑ Dispense/CLICK: 0.25 mL
- ☑ Approx. Residual: 2.5 mL
- ☑ Approx. Priming: 0.25 mL or 1 CLICK

USA | CANADA | AUSTRALIA PACIFIC | EUROPE
WWW.MEDISCA.COM